

Dr's Cam



-IC-DCAM70 -IC-DCAU71

Enjoy Wireless & Auto Focus!

Wireless Autofocus Intraoral Camera with Mouse

1 Features of Dr's cam

Wireless type only



High speed digital wireless

High speed digital wireless provide more clear live image without interruption.



Wireless receiver

Compact wireless receiver without any power supply.



Auto focus

Can get clear image with auto focus function in oral cavity.



HD image

Can get high definition image without distortion.



Gyro mouse built-in

eable as a mouse.



280° twist head

Easy to get approach anywhere in oral cavity.



OLED screen

Bright and clear OLED screen indicates various device status.



Adjustable LED brightness

To avoid image quality loss by light reflection, can adjust LED brightness.



3 Color Impressions

Can change live image color impression (Original, Warm, Cool)



Software compatibility

Can be use camera with most of dental imaging software. * It would need twain driver for some of dental imaging software.



2 Specification

Model name		IC-DCAM70(Wireless model)
Camera	Operation voltage	3.7V
	Size	Approx. 255 x 30 x 28 mm
	Weight	Approx. 89g
	Resolution	1.3 Mega pixel
	Sensor/Lens type	CMOS / Liquid lens
Charger	Operation voltage	5V
	Size	Approx. 110 x 105 x 68 mm
	Weight	Approx. 100g
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

Model name		IC-DCAU/1(Wired model)
Camera	Operation voltage	5V
	Size	Approx. 202 x 30 x 28 mm
	Weight	Approx. 56g
	Resolution	1.3 Mega pixel
	Sensor/Lens type	CMOS / Liquid lens
USB cable	Length	Approx. 3M
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

3 Minimum System Requirement

PC system:	2.0 Ghz CPU, 1GB RAM	
os:	Windows XP service pack 3	
USB port:	USB 2.0	









Upper images are take by real product.

4 DrsViewer2

Image Viewer included.







Patient Registration Live View

Image Edit

Print View





Good\(\bar{\D}\)rs

Tel: +82 32 424 6325 Fax: +82 32 424 6326 E-mail: good@gooddrs.com Website: www.gooddrs.com